


JET ENGINE TRAINER



Experimental capabilities

- Identification of the components
- Standard Operation Procedure of a Jet Engine
- Analysis of the operation of a jet engine
- Determination of the specific thrust
- Determination of the specific fuel consumption
- Determination of turbine efficiency
- Temperature measurement at different points
- Pressure measurement at different points
- Security system for the operation of the turbine
- Recording of characteristics curves
- Diagramm T-S (temperature in function of entropy)

Operating principle

The TGT 020 trainer investigates the operation of a jet engine.

TGT 020 includes the following components: jet engine (with compressor, combustion chamber, turbine and propelling nozzle), fuel system, starter and ignition system and measuring and control equipment. The gas turbine works as an open cyclic process, with the ambient air being drawn out and fed back in the jet engine the ambient air drawn in is first brought to a higher pressure in the single-stage radial compressor. In the downstream combustion chamber, fuel is added to the compressed air and the resulting mixture is combusted. The temperature and flow velocity of the gas increase. The gas flows out of the combustion chamber into the single-stage axial turbine and discharges a portion of its energy to the turbine. This turbine drives the compressor. In the propelling nozzle, the partially expanded and cooled gas expands to ambient atmospheric pressure and the gas accelerates to almost the speed of sound. The high-speed gas outflow generates the thrust. In order to reduce the outlet temperature, the exhaust gas stream is mixed with the ambient air in a mixing pipe.

Between the compressor and the turbine is the annular combustion chamber. With optimum fuel utilization, low pressure loss and good ignition properties, the ring shape of this combustion chamber is typical for the design used in jet engines. The movable turbine platform, with a strain gauge, enables measurement of the thrust.

The speed, temperatures and mass flow rates of the air and fuel are recorded using sensors.


The measured values are transmitted directly to a computer. The data acquisition software is included.

The robust design of this device makes it suitable for use in schools.

The equipment is set up on an aluminium profile. The frame is also equipped with adjustable feet with rubber to avoid vibrations while the engine is running. This gives it great strength and a flexibility of integration into your laboratory.

The manufacture of this equipment complies with the European standard for machinery manufacturing.

Illustrations



1) Fuel tank

The turbine can use the following fuel :

- kerosene JET-A1 + turbine oil
- petroleum de-aromatized + turbine oil

Volume : 20L

Drain valve at the bottom and vent and filling hole on the top.

The tank is made of stainless steel and is equipped with a lateral level gauge

2)Electrical box

The machine includes an electrical box conforming to European standards. It contains at least:

- a general power disconnecter
- a 30mA differential circuit breaker
- relaying and circuit breakers necessary for operation
- the button and indicators necessary for operation
- an emergency stop button
- a 7" touch screen to display the measure locally
- a potentiometer to adjust the power of the jet engine
- the pocket controller for the jet engine

3) Frame

Transparent protective screen for better visibility and safety

The unit is set up on an anodized aluminum frame on castors wheels

Technical details

4) Jet engine

Brand : Jet Cat

Type : P200RX

Max thrust : 230 N at 112 000 min⁻¹

Speed range : 33 000 to 112 000 min⁻¹

Fuel consumption : max 850 ml/min

Exhaust gas temperature : 200 °C at 1 m distance

Sound level at 1 m distance : 130 dB

Electrical start of turbine

5)Measurements

The system is fully instrumented to collect the data and study the jet engine. All the data are displayed on a local touch screen and are sent to the data acquisition system as well.

The system includes the following sensors :

Pressure sensors (digital sensors):

atmospheric pressure (P_{atm}), air inlet (P₁), outlet of the compressor (P₂), inlet of the turbine (P₃), outlet of the turbine (P₄)

Temperature sensors (thermocouple K type):


atmospheric temperature (T_{atm}), fresh air inlet (T₁), outlet of the compressor (T₂), fuselage (T₃), inlet of the turbine (T₄), outlet of the turbine (T₅), nozzle outlet (T₆)

Air velocity (pitot tube) : (m/s)

Rotational Speed of the turbine (optical system): (rpm)

Thrust (load cell) : (kg)

Fuel consumption (turbine flowmeter): (cm³/min)



TGT020



DATA ACQUISITION SYSTEM

The trainer includes as standard data acquisition software and control system. The connection towards the PC is made via a standard Ethernet port. All these values can be saved to an Excel file to allow further analysis. This application can be installed on any computer (latest version of windows) and is free of license


IDENTIFICATION :

There is in this window the illustration of the turbine to identify the different part of turbojet



SYNOPTIC:

There is in this window the block diagram of the jet engine with the location of the measurement points and their values.



GRAPH :

There is in this graph window, the possibility of drawing the measurement curves as a function of the time by selecting the desired measure.





DIAGRAM T-S :


We find in this graph window, the possibility to draw a curve with two parameters. For example, the temperature depending on the entropy. The curve is drawn step by step depending on the evolution of the process.



CONTROL OF THE ENGINE :

The software can also be used to control the jet engine directly from the computer. Those actions can be done:

- The control of the turbine (start, stop..)
- Control of the power (0-100%)



Services required

- Electrical supply: 230 VAC - 50 Hz – 20 A
- Fuel : kerosene Jet-A1, or petroleum dearomatised + turbine oil mix at 5% (1 L oil included)
- Exhaust gases extractor in stainless steel
- Dimensions: (L x W x H mm): 1100 x 785 x 1780
- weight (Kg): 130

Note : if the equipment installation is operated by our staff, all supplies and exhaust connections required must stand at less than 2m from the machine

Documentation

- User's manual
- Technical documentation of the components
- Lab exercises
- Wiring diagram
- Data acquisition software
- Certificate of conformity CE

DIDATEC– Zone d'activité du parc – 42490 FRAISSES- FRANCE
Tél. +33(0)4.77.10.10.10 – Fax+33(0)4.77.61.56.49 – www.didatec-technologie.com
email : service_commercial@didatec-technologie.com

Reproduction interdite / copy prohibited– Copyright DIDATEC mai-20- page 3

Dans le cadre de l'amélioration permanente de nos produits, ce descriptif technique est susceptible d'être modifié sans préavis
As part of the continuous improvement of our products, this technical specification may be modified without previous notifying

Illustrations non contractuelles / Illustrations not contractual

version : FT-TGT020-STD-C