CRANH3



REFRIGERATION CYCLE BY ABSORPTION



Experimental capabilities

- Study and basic concept of a refrigeration installation by absorption of scattering type.
- Study of the important phases of the cycle.
- Analysis of the heat balances in the evaporator and in the other components.

CRANH3



Operating principle

An absorption refrigeration system uses three products: hydrogen, water and ammonia. A small amount of sodium chromate Na2CrO4 is added to reduce corrosion, but this product has no effect on cooling.

A solution of ammonia (ammonia in solution in water) is boiled in a boiler, the temperature is about 180 °. Ammonia that boils at a lower temperature than water escapes from the solution. Boiling carries away part of the solution (which is here a very dilute solution because practically all the ammonia has vaporized).

The water is harvested in a separator and returned to the boiler. The temperature is about 110° and the pressure causes the remaining water vapour to liquefy. The boiler has a peculiar shape: a boiling tube in the center and a siphon on the outside. The water that falls back into the boiler actually falls into the outer siphon. Boiling ammonia carries water droplets to the siphon (this is the 'bubble pump').

The ammonia that evaporates forms bubbles that escape from the solution and carry ammonia (more and more diluted). We have ammonia gas and ammonia water at the outlet of the bubble pump.

The ammoniacal water flows into the siphon and then into the absorber via a hose thanks to gravity. In the pipe, the water is cold The ammonia cools and liquefies in the condenser. Ammonia comes into contact with hydrogen in the evaporator. The partial pressure of ammonia decreases and ammonia evaporates, producing cold.

The gaseous ammonia now passes through a tank and rises to the absorber where it encounters water that runs against the current. The ammonia gas dissolves in the water and we again obtain a concentrated solution of ammonia. As a gas, only hydrogen remains, which is again used in the evaporator. The absorption of ammonia gas in runoff water keeps the partial pressure of ammonia low.

The rugged design of this equipment makes it perfectly suited for use in schools.

Its anodized aluminum structure on wheels gives it great robustness and flexibility of integration into your premises. The manufacture of this equipment complies with the European machinery directive.

Illustrations

Fusible de securité Echangeur de température du gaz Absorbeur Absorbeur Récipient de l'absorbeur de l'absorbeur du liquide Solution faible en ammoniaque Solution forte en ammoniaque Uquide d'ammoniaque Gaz hydrogène Wapdragène Wydrogène et vapeur d'ammoniaque Gaz hydrogène Wydrogène Wydrogè

Technical details

- The equipment consists of a generator, an absorber, a condenser, and a tank;
- An instrumentation console will make it possible to raise the power consumed over time, the load as well than temperatures;
- 3. Temperature indicators at different characteristic points (8)
- 4. Power indicator 0-150W
- 5. Variable electrical load
 - Electric heater: 125 W
 - Variable heating generator
 Evaporator heating, adjustable: 50 W
 - Adjustable gas burner: propane
- 7. The unit is protected by switches, circuit breakers as well as a differential circuit breaker to protect the user
- 8. Working fluid: ammonia-water
- 9. Auxiliary gas: hydrogen

Services required

- Electrical supply: 230 Vac 50 Hz
- Propane or butane gas supply (cylinder)
- Dimensions: (LxWxH mm): 700 x 600 x 1100
- weight (Kg): 60

Note: if the equipment installation is operated by our staff, all supplies and exhaust connections required must stand at less than 2m from the machine

Documentation

- Instructions
- Pedagogical manual
- Technical file
- Practical work
- Electrical diagram
- Fluidic diagram
- · EC certificate of conformity

DIDATEC- Zone d'activité du parc - 42490 FRAISSES- FRANCE Tél. +33(0)4.77.10.10.10 - Fax+33(0)4.77.61.56.49 - <u>www.didatec-technologie.com</u>